

QA3™

Treatment Stretcher

Supreme comfort, optimal positioning and optimal access for a wide range of clinical treatments – the QA3 Treatment Stretcher offers both the user and the patient an improved clinical journey

Powered ⚡ Functions

Powered ⚡
Backrest

Powered ⚡
Trendelenburg

Powered ⚡
Height adjustment

Hand Control +
Foot-end Control

Powered ⚡
Profile (Back & Knee)

Powered ⚡
KneeFlex

Powered ⚡
Reverse Trendelenburg

Powered ⚡ CPR



For more information on any of our products or service contracts, call:

sales@aneticaid.com

+44 (0)1943 878647

aneticaid.com

SPECIFICATIONS			
Minimum height	*495 mm	Mattress length	2025 mm
Maximum height	*875 mm	Mattress width	700 mm
*excludes mattress		Mattress depth	90 mm
Minimum length	1700 mm	Maximum patient weight	320kg
Maximum length	2150 mm	Backrest	0 to +90°
Overall width	905 mm	KneeFlex	0 to +30°
Weight	141 kg	Trendelenburg	12° to -12°

Innovative **Medical** Technology - **Practically Applied**

Anetic Aid Ltd, Sapper Jordan Rossi Business Park, Otley Road, Baildon, BD17 7AX, UK

AneticAid